

# Beef Clip Rodder



## Beef Rodder

The Adept Beef Rodder is used to separate the weasand (oesophagus) away from the trachea and from other connecting tissue. The Adept Beef Weasand Clip is then applied and pushed up into position using the Adept Beef Rodder.

- Designed for optimum separation of connecting tissue without cutting
- Lightweight and ergonomic
- Polished finish for ease of cleaning and sterilization
- Fast and easy to use
- Available in two sizes

Product Code	0142-010 – Beef Clip Rodder (White) 0142-030 – Beef Clip Rodder (Orange) 0142-040 – Beef Clip Rodder (Green) 0142-070 – Beef Clip Rodder (Long)
Colours	White, Orange, Green
Material	Stainless Steel Plastic Handle
Dimensions	0142-010 - 675mm 0142-030 - 675mm 0142-040 - 675mm 0142-070 - 825mm